



EPI FIX INSTRUCTION LEAFLET



CE



CHRIS HANLEY & PARTNERS
4 Kilvey Road, Brackmills, Northampton NN4 7BQ
T: 01604 709999 F: 01604 700057
www.chaneco.co.uk



Features

- Moulded and shaped thermoplastic clasp for comfort and versatility
- Padded band places pressure directly onto the inflamed tendon
- Silicone dots provide improved suspension onto the skin
- Easily adjusted by elastic band closure
- Promotes healing and reduces pain
- Low profile, lightweight and durable

Indications

- Tennis Elbow (Lateral Epicondylitis)
- Golfers Elbow (Medial Epicondylitis)
- Prevention of sporting / work injuries
- Treatment of repetitive strain injuries
- Bursitis

Patient Information

- This product should be used only on the recommendation of a suitably qualified medical professional
- The initial fitting should be carried out by a suitably qualified medical professional
- You should follow the instructions and advice given to you regarding use of this product

Fitting Instructions

1. Place around the arm with plastic loop on the top of your arm, 3cm (approximately 2 fingers) down from the elbow to allow complete bending to 90°
2. Ensure that the pad is in the correct position to exert pressure on the tendon
3. Adjust the Velcro strap to achieve the desired traction using the red dots as a guide. Do not exceed the maximum limit (solid red section)
4. Add or remove the extra white spacer to achieve greater or lesser compression as required

Care and Maintenance

- Stick the Velcro pads together before washing
- Hand wash only in warm water and mild detergent
- Use a dry towel to absorb excess moisture and then allow to dry naturally, away from direct heat sources. Do not iron
- Do not use alcohol, creams or liquid solvents during use or cleaning

Single Patient Use

Any product supplied by Chris Hanley & Partners is for use by a single named patient only, as prescribed and fitted by a qualified clinician. Re-use of any product is potentially dangerous and will invalidate the CE mark and any warranty.