

CE



CHRIS HANLEY & PARTNERS
4 Kilvey Road, Brackmills, Northampton NN4 7BQ
T: 01604 709999 F: 01604 700057
www.chaneco.co.uk



FIXED SWEDISH KNEE CAGE INSTRUCTION LEAFLET



Introduction

An effective Orthosis used in the management of Hyperextension, the Swedish Knee Cage is quick and easy to fit. An adjustable posterior pressure pad, attached to the upright, allows for an individual fit, whilst the straps are lined with a strong elastic for a firm grip of the leg.

- Provides management of genu recurvatum
- Lightweight vinyl coated aluminium moulds easily
- Adjustable posterior pressure pad
- Permits flexion at knee



Fitting Instructions

1. Unfasten one side of the proximal and distal straps at the snaps and slide off the uprights
2. The patient should be seated with the knee in flexion
3. Apply the brace from the posterior, centering the middle back pad to the knee
4. Adjust the side sections of the brace to the knee for comfort
5. Bring the proximal and distal straps across the front of the leg and reattach to the uprights. The elastic section of the straps should be against the leg.
6. Check the degree of knee extension with the patient standing. Adjust if necessary.

NOTE: The position of the back pad controls the amount of extension lock-out i.e. A more anterior position provides less lock-out (more knee extension)

Single Patient Use

Any product supplied by Chris Hanley & Partners is for use by a single named patient only, as prescribed and fitted by a qualified clinician.

Re-use of any product is potentially dangerous and will invalidate the CE mark and any warranty.