affix

REF.20050

Peroneusorthese
Ankle foot orthosis
Tuteur pour pied ballant
Тутор на голеностопный сустав
Tutore per piede ciondolante

ORTHOSERVICE
RO+TEN

Headquarter: ORTHOSERVICE AG
Via Milano 7 · CH-6830 Chiasso (TI) · Switzerland
Tel. 0041 (0) 91 822 00 88 · Fax 0041 (0) 91 822 00 89
info@orthoservice.com · www.orthoservice.com
Niederlassung Deutschland: Orthoservice Deutschland GmbH
Flugstraße 8 · D-76532 Baden-Baden · Deutschland
Tel. 0049 (0) 7221 991 39 11 · Fax 0049 (0) 7221 991 39 13
info@orthoservice.de · www.orthoservice.de
Sede italiana: RO+TEN s.r.l.
Sede legale: Via Marco De Marchi, 7 · I-20121 Milano (MI) · Italia
Sede operativa e amministrativa:
Via Comasina, 111 · I-20843 Verano Brianza (MB) · Italia
Tel. 0039 039 601 40 94 · Fax 0039 039 601 42 34
info@roplusten.com · www.roplusten.com
Società soggetta a Direzione e Coordinamento (art. 2497bis CC):
Orthoservice AG (CH) · 6830 Chiasso (TI) · Switzerland

DEUTSCH

ENGLISH

FRANÇAIS

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ITALIANO

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Please read the instructions carefully and thoroughly

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares with sole responsibility that this product is a Class I medical product and has been manufactured in accordance with the basic requirements of the directives 93/42/EEC and 2007/47/EEC. These operating instructions were created in accordance with the specification of the above-mentioned directives. They aim to ensure the safe and correct use of the medical product.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SELECTION/SIZE

Code: REF. 20050				
Size	XS	S	M	L
Ankle circumference cm	13/17	17/21	21/25	25/29
Colour: black				
Fits right and left				

SELECTION/SIZE ACCESSORY

Code: REF. 20051			
Size	S	M	L
Midfoot circumference cm	17/21	21/25	25/29
Colour: black			
Specify right or left			

MAINTENANCE

- Do not bleach
- Do not iron
- No chemical cleaning
- Do not tumble-dry
- Washing instructions: hand wash at max. 30° C with neutral soap. Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

PRECAUTIONS

It is advised that the pressure exerted by the product does not act on parts of the body suffering from wounds, swelling or lumps. The product should not be tightened too firmly so as to avoid the development of pressure points or compression of the nerves or blood vessels underneath. It is advisable to wear the product over a garment, avoiding direct contact with the skin. If in doubt about how to apply, contact an orthopaedic technician. The device should not be worn in the vicinity of open flames or strong electromagnetic fields.

WARNINGS

The product may only be used in the event of medical indication or when prescribed by a physician and must be applied by an orthopaedic technician, the competent authority for the application as well as for information on the safe use, according to individual needs. To ensure the effectiveness, tolerability, and correct operation, it is necessary that the device is put on with the utmost care. Any changes to the structure or adjustment of the device must be prescribed by a physician and performed by an orthopaedic technician. Under no circumstances should the adjustment made by the physician/orthopaedic technician be altered. Use of the product is recommended for a single patient. With hypersensitive patients, direct contact with skin may cause redness or irritation. In the case of onset of pain, swelling, lumps or any unusual reaction, immediately contact your physician.

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ENGLISH

CHARACTERISTICS AND MATERIALS

- To be worn ONLY with lace-up shoes.
- Low profile.
- Breathable micropore fabric.
- A soft pad attaches the support to the shoe.
- Tensioning elastic regulates the degree of flexion.
- Gel padding over the Achilles tendon.
- Fastens with Velcro® straps.
- Additional midfoot band for wearing with bare feet (optional).

INDICATIONS

- Flaccid paralysis of the foot
- Light to moderate spastic equinus foot

CONTRAINDICATIONS

Currently no known

PUTTING ON THE APPLIANCE

- A** Separate the three parts of the support (shoe pad, tensioning elastic and the tibial cuff) and fully open the tibial cuff.
- B** With the tibial cuff open, rest it above the back of the ankle, with the gel padding over the Achilles tendon and the bend upwards (fig.1a), and wrap it around the ankle using the two Velcro® straps (fig.1b). Leave the D-ring completely free and uncovered.
- C** Insert the shoe pad between the laces and the tongue of the shoe (fig.2).
- D** Put on the shoe and lace it up (fig.3).
- E** With the tensioning elastic folded in two, pass it through the D-ring from top to bottom (fig.4).
- F** Attach the folded end of the elastic to the shoe pad (fig.5).
- G** Finally, attach the elastic to the support and achieve the desired degree of tension by fastening the two Velcro® ends level with the malleoli (fig.6).

ACCESSORY: MIDFOOT BAND

CHARACTERISTICS AND MATERIALS

- To be worn with bare feet or open shoes together with the Affix.
- Low profile.
- Soft, breathable micropore fabric.
- Anatomically designed to fit the left or right foot.
- Fastens with Velcro®.

PUTTING ON THE APPLIANCE

- A** With the tibial cuff open, rest it above the back of the ankle, with the gel padding over the Achilles tendon and the bend upwards (fig.10a), and wrap it around the ankle using the two Velcro® straps (fig.10b). Leave the D-ring completely free and uncovered.
- B** Fit the band onto the midfoot with the ring on the top (fig.11).
- C** With the tensioning elastic folded in two, pass it through the D-ring from top to bottom (fig.12).
- D** Attach the folded end of the elastic to the ring on the midfoot band (fig.13).
- E** Finally, attach the elastic to the support and achieve the desired degree of tension by fastening the two Velcro® ends level with the malleoli (fig.14).

WARRANTY

This warranty is governed by the statutory provisions of the country in which the product was purchased. Where relevant that country's regulations concerning warranties and the relationship between the seller and buyer shall apply. If you believe that the terms of the warranty may be applicable, first contact the retailer where the product was purchased. Orthoservice will repair or replace all or part of the device and its accessories if defects in materials and workmanship arise during the period of six months from the date of purchase.

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