

Care and Maintenance

- Wash periodically by hand with a damp cloth and mild soap
- Use a towel to absorb most of the dampness and then allow to dry at room temperature.
- Do not hang up, or iron and do not expose to direct heat sources such as stoves, heaters, radiators, direct sunlight etc
- Do not expose to alcohol, ointments or solvents.



Achilles Wedge Fitting Guidelines & Care Instructions



CS Lifesciences Europe Limited,
The Black Church,
St Mary's Place, Dublin 7, Dublin,
D07 P4AX, Ireland

Version 4
Nov 2021



CHANECO LTD
4 Kilvey Road, Brackmills, Northampton, NN4 7BQ
T: 01604 709999 F: 01604 700057
www.chaneco.co.uk



Single Patient Use



Not Fire
Retardant



Chaneco Limited is registered in England under Company Number 11511502. Registered office as above.

Features

- 4" (10cm) neoprene padded foam wedge
- Fits into most closed heel walkers
- Multi-layered for progressive tendon stretching
- Each section is $\frac{3}{4}$ " (2cm) in height and may be "peeled away" during Achilles Tendon rehabilitation to achieve the desired height.

Indications

- Achilles Tendon Rehabilitation
- Post-Trauma
- Post-Operative

Size Table

Size	Code
Universal	AW-1

Fitting Instructions

- Peel away sections of the wedge to reach the desired angle as required
- Place wedge into the walker
- Remove another section every two weeks while rehabilitation progresses.
- Patient can then use the single, smallest wedge in their own footwear for one or two weeks
- Always following the guidance from your physician or Orthotist when using this product.

Additional Information

- Contains Latex