hallus valgus

REF. 360

Hallux-valgus-bandage
Hallux valgus retractor
Ecarteur pour hallux valgus
Орtez для большого пальца
стопы при бурсите
Divaricatore per alluce valgo

DEUTSCH

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ENGLISH

Please read the instructions carefully and thoroughly

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG states under its sole responsibility that this product is a Class I medical device that has been manufactured in accordance with the requirements of Directives 93/42/EEC and 2007/47/EC. These instructions have been written in conformity with these directives, and their aim is to ensure the safe and proper use of this medical device.

SELECTION/SIZE

Size Shoe size Colour	S 34/36 white	M 37/40	L 41/46
specify right or left			

MAINTENANCE

- Do not bleach No chemical cleaning
- Do not iron Do not tumble-dry
- Washing instructions: hand wash at max. 30° C with neutral soap.
Dry away from heat sources. Clean with a sponge soaked in handwarm water (30° C) and neutral soap. Dry with a cloth.

Do not dispose of the product or any of its components into the environment

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

PRECAUTIONS

The pressure exerted by the device should not act on the parts of the body which have injuries and swelling. A sock should be worn. Highly sensitive patients may suffer from a reddening of the skin and irritation if there is direct contact.

WARNINGS

The product may only be used in the event of medical indication or when prescribed by a physician and must be applied by an orthopaedic technician, the competent authority figure for the application as well as for information on the safe use, according to individual needs.
To ensure the effectiveness, tolerability, and correct operation, it is necessary that the device is put on with the utmost care.
Any changes to the structure or adjustment of the device must be prescribed by a physician and performed by an orthopaedic technician. Under no circumstances should the adjustment made by the physician/orthopaedic technician be altered.
Use of the product is recommended for a single patient.
With hypersensitive patients, direct contact with skin may cause redness or irritation. In the case of onset of pain, swelling, lumps or any unusual reaction, immediately contact your physician.

WARRANTY

This warranty is governed by the statutory provisions of the country in which the product was purchased. Where relevant that country's regulations concerning warranties and the relationship between the seller and buyer shall apply. If you believe that the terms of the warranty may be applicable, first contact the retailer where the product was purchased. Orthoservice will repair or replace all or part of the device and its accessories if defects in materials and workmanship arise during the period of six months from the date of purchase.

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