

Care and Maintenance

- Wash periodically by hand with a damp cloth and mild soap
- Use a towel to absorb most of the dampness and then allow to dry at room temperature.
- Do not hang up, or iron and do not expose to direct heat sources such as stoves, heaters, radiators, direct sunlight etc
- Do not expose to alcohol, ointments or solvents.



Hinged Knee Wrap Fitting Guidelines & Care Instructions



KSH IFU - Revision 4 - Jan 2024

CHANECO LTD
4 Kilvey Road, Brackmills, Northampton, NN4 7BQ
T: 01604 709999 F: 01604 700057
www.chaneco.co.uk

Chaneco Limited is registered in England under Company Number 11511502. Registered office as above.



Single Patient Use



Not Fire
Retardant



Features

- Front opening, wrap-around design with enclosed hinged side steels for increased medial/lateral stability
- Provides compression and therapeutic warmth while the open patella and popliteal area add comfort during flexion
- Hinged side steels easily removed through open pockets
- Secure hook and loop fastening
- Available in either black Ortho-Tex or grey Neoprene and in two lengths - Standard or Long (Long length available in Neoprene only)
- Standard length 12" (30cm), Long length 15¾" (40cm)
- Universal - fits left and right knee

Indications

- R/A and O/A
- Mild to moderate hyperextension
- Knee instability

Size Table

Size	Knee Circumference		Standard		Long Neoprene
	Inches	CM	Ortho-Tex	Neoprene	
X Small	12½ - 13½	32 - 34	KSH-S-O-XS	KSH-S-N-XS	KSH-L-XS
Small	13½ - 14½	34 - 37	KSH-S-O-S	KSH-S-N-S	KSH-L-S
Medium	14½ - 15	37 - 38	KSH-S-O-M	KSH-S-N-M	KSH-L-M
Large	15 - 16	38 - 41	KSH-S-O-L	KSH-S-N-L	KSH-L-L
X Large	16 - 17	41 - 43	KSH-S-O-XL	KSH-S-N-XL	KSH-L-XL

Fitting Instructions

1. Undo the tabbed Velcro straps and open up the wrap.
2. Wrap around the leg with the inside label at the top and the popliteal opening at the back, making sure that the hinges line up with the knee cap and that it is visible through the opening
3. Fasten the wrap using the Velcro straps to ensure a snug fit around the leg.

Additional Information

- Do not use if heat causes allergic reactions in the patient
- LATEX FREE