

## Care and Maintenance

- Wash periodically by hand with a damp cloth and mild soap
- Use a towel to absorb most of the dampness and then allow to dry at room temperature
- Do not dry clean or tumble dry - allow to dry naturally away from direct heat sources
- Do not expose to alcohol, ointments or solvents



## Triple Strap Maternity Belt

### Fitting Guidelines & Care Instructions



MBTS IFU - Revision 4 - Jan 2024

**CHANECO LTD**  
4 Kilvey Road, Brackmills, Northampton, NN4 7BQ  
T: 01604 709999 F: 01604 700057  
[www.chaneco.co.uk](http://www.chaneco.co.uk)

Chaneco Limited is registered in England under Company Number 11511502. Registered office as above.



Single Patient Use



Not Fire Retard-



## Features

- Elastic 3-strap design with secure hook and loop closure
- Allows the patient to obtain optimum fit through a range of strap positions, supporting the abdomen and reducing strain on the pelvis
- Strong lightweight elastic gives excellent support and comfort without being too restrictive
- Expands with increase in abdominal growth
- Easy to put on, take off and adjust whilst wearing

## Indications

- Pubic Symphysis pain
- Sacro-iliac pain
- Helps relieve lower back pain during pregnancy
- Can be used during post-natal recovery

## Size Table

Size	Abdomen Circ.		Code
	Inches	CM	
Small	35½ - 40¼	90 - 102	MBTS-S
Medium	40¼ - 45¾	102 - 116	MBTS-M
Large	45¾ - 53¼	116 - 136	MBTS-L
X Large	53¼ - 61¼	136 - 156	MBTS-XL

## Fitting Instructions

- Place the belt over the sacro-iliac joints at the back and wrap the straps around the pelvis, directing them forwards and downwards to lie under the lower abdomen and across the pubic bone.
- You can position the three straps all together on top of each other for greater support or place the three strap in different positions to cover a wider range of support if required.
- Close the belt at the front using the hook and loop fastenings.

## Additional Information

- Material content: 70% Elastic strap, 20% cotton, 10% nylon
- NOT LATEX FREE