

Ladies - CITY, CONTEMPORARY & CLASSIC RANGE - LINEA LAST - Fitting Chart C

City, Contemporary & Classic Range Linea Lasts have a 3mm Inlay Allowance as Standard
 Footwear is supplied with 1 x 3mm Poron Inlay & 1 x Polyurethane Footbed as Standard
 (For 1/2 Sizes use measures half way between whole sizes)

Ladies Slim Fitting

Regular Lasts						Alice Lasts (Narrow Heel)		Chloe Lasts (Shallow Forepart)		Emma Lasts (Narrow Heel & Shallow Forepart)			
Size	Heel Width (mm)	Joint Width (mm)	Joint Circ. (mm)	Instep Circ. (mm)	Toe Depth (mm)	Heel Width (mm)	Instep Circ. (mm)	Joint Circ. (mm)	Toe Depth (mm)	Heel Width (mm)	Joint Circ. (mm)	Instep Circ. (mm)	Toe Depth (mm)
1	53	79	192 - 204	180 - 192	21 - 27	49	172 - 184	184 - 196	17 - 23	49	184 - 196	172 - 184	17 - 23
2	54	82	199 - 211	187 - 199	22 - 28	50	179 - 191	191 - 203	18 - 24	50	191 - 203	179 - 191	18 - 24
3	56	84	206 - 218	193 - 205	23 - 29	52	185 - 197	197 - 209	18 - 24	52	197 - 209	185 - 197	18 - 24
4	57	87	212 - 224	200 - 212	23 - 29	53	192 - 204	204 - 216	19 - 25	53	204 - 216	192 - 204	19 - 25
5	58	89	218 - 228	206 - 212	24 - 30	55	199 - 211	210 - 222	20 - 26	55	210 - 222	199 - 211	20 - 26
6	60	92	225 - 237	213 - 225	25 - 31	56	205 - 217	217 - 229	21 - 27	56	217 - 229	205 - 217	21 - 27
7	61	94	231 - 243	219 - 231	26 - 32	58	212 - 224	223 - 235	22 - 28	58	223 - 235	212 - 224	22 - 28
8	63	97	238 - 250	226 - 238	26 - 32	59	218 - 230	230 - 242	23 - 29	59	230 - 242	218 - 230	23 - 29
9	64	99	244 - 256	232 - 244	27 - 33	61	225 - 237	236 - 248	23 - 29	61	236 - 248	225 - 237	23 - 29
10	66	102	251 - 263	239 - 251	28 - 34	62	231 - 243	243 - 255	24 - 30	62	243 - 255	231 - 243	24 - 30

(Any measures not given for our Alice, Chloe or Emma Lasts are as our Regular Lasts)

Ladies Standard Fitting

Regular Lasts						Alice Lasts (Narrow Heel)		Chloe Lasts (Shallow Forepart)		Emma Lasts (Narrow Heel & Shallow Forepart)			
Size	Heel Width (mm)	Joint Width (mm)	Joint Circ. (mm)	Instep Circ. (mm)	Toe Depth (mm)	Heel Width (mm)	Instep Circ. (mm)	Joint Circ. (mm)	Toe Depth (mm)	Heel Width (mm)	Joint Circ. (mm)	Instep Circ. (mm)	Toe Depth (mm)
1	57	85	204 - 216	188 - 200	21 - 27	53	180 - 192	196 - 208	17 - 23	53	196 - 208	180 - 192	17 - 23
2	58	88	210 - 222	195 - 207	22 - 28	54	187 - 199	203 - 215	18 - 24	54	203 - 215	187 - 199	18 - 24
3	60	90	217 - 229	201 - 213	23 - 29	56	193 - 205	209 - 221	18 - 24	56	209 - 221	193 - 205	18 - 24
4	61	93	223 - 235	208 - 220	23 - 29	57	200 - 212	216 - 228	19 - 25	57	216 - 228	200 - 212	19 - 25
5	63	95	230 - 242	214 - 226	24 - 30	59	206 - 212	222 - 234	20 - 26	59	222 - 234	206 - 212	20 - 26
6	64	98	236 - 248	221 - 233	25 - 31	60	213 - 225	229 - 241	21 - 27	60	229 - 241	213 - 225	21 - 27
7	66	100	243 - 255	227 - 239	26 - 32	62	219 - 231	235 - 247	22 - 28	62	235 - 247	219 - 231	22 - 28
8	67	103	250 - 262	234 - 246	26 - 32	63	226 - 238	242 - 254	23 - 29	63	242 - 254	226 - 238	23 - 29
9	69	105	256 - 268	240 - 252	27 - 33	65	232 - 244	248 - 260	23 - 29	65	248 - 260	232 - 244	23 - 29
10	70	108	263 - 275	247 - 259	28 - 34	66	239 - 251	255 - 267	24 - 30	66	255 - 267	239 - 251	24 - 30

(Any measures not given for our Alice, Chloe or Emma Lasts are as our Regular Lasts)

Ladies Wide Fitting

Regular Lasts						Alice Lasts (Narrow Heel)		Chloe Lasts (Shallow Forepart)		Emma Lasts (Narrow Heel & Shallow Forepart)			
Size	Heel Width (mm)	Joint Width (mm)	Joint Circ. (mm)	Instep Circ. (mm)	Toe Depth (mm)	Heel Width (mm)	Instep Circ. (mm)	Joint Circ. (mm)	Toe Depth (mm)	Heel Width (mm)	Joint Circ. (mm)	Instep Circ. (mm)	Toe Depth (mm)
1	61	91	216 - 228	196 - 208	21 - 27	57	188 - 200	208 - 220	17 - 23	57	208 - 220	188 - 200	17 - 23
2	62	94	222 - 234	203 - 215	22 - 28	58	195 - 207	214 - 226	18 - 24	58	214 - 226	195 - 207	18 - 24
3	64	96	229 - 241	209 - 221	23 - 29	60	201 - 213	221 - 233	18 - 24	60	221 - 233	201 - 213	18 - 24
4	65	99	235 - 247	216 - 228	23 - 29	61	208 - 220	227 - 239	19 - 25	61	227 - 239	208 - 220	19 - 25
5	67	101	242 - 254	222 - 234	24 - 30	63	214 - 226	234 - 246	20 - 26	63	234 - 246	214 - 226	20 - 26
6	68	104	248 - 260	229 - 241	25 - 31	64	221 - 233	240 - 252	21 - 27	64	240 - 252	221 - 233	21 - 27
7	70	106	255 - 267	235 - 247	26 - 32	66	227 - 239	247 - 259	22 - 28	66	247 - 259	227 - 239	22 - 28
8	71	109	262 - 274	242 - 254	26 - 32	67	234 - 246	254 - 266	23 - 29	67	254 - 266	234 - 246	23 - 29
9	73	111	268 - 280	248 - 260	27 - 33	69	240 - 252	260 - 272	23 - 29	69	260 - 272	240 - 252	23 - 29
10	74	114	275 - 287	255 - 267	28 - 34	70	247 - 259	267 - 279	24 - 30	70	267 - 279	247 - 259	24 - 30

(Any measures not given for our Alice, Chloe or Emma Lasts are as our Regular Lasts)

Ladies Deep Fitting

Regular Lasts					
Size	Heel Width (mm)	Joint Width (mm)	Joint Circ. (mm)	Instep Circ. (mm)	Toe Depth (mm)
1	57	85	216 - 228	188 - 200	27 - 33
2	58	88	222 - 234	195 - 207	28 - 34
3	60	90	229 - 241	201 - 213	29 - 35
4	61	93	235 - 247	208 - 220	29 - 35
5	63	95	242 - 254	214 - 226	30 - 36
6	64	98	248 - 260	221 - 233	31 - 37
7	66	100	255 - 267	227 - 239	32 - 38
8	67	103	262 - 274	234 - 246	32 - 38
9	69	105	268 - 280	240 - 252	33 - 39
10	70	108	275 - 287	247 - 259	34 - 40

Ladies Wide & Deep Fitting

Regular Lasts					
Size	Heel Width (mm)	Joint Width (mm)	Joint Circ. (mm)	Instep Circ. (mm)	Toe Depth (mm)
1	61	91	228 - 240	196 - 208	27 - 33
2	62	94	235 - 247	203 - 215	28 - 34
3	64	96	241 - 253	209 - 221	29 - 35
4	65	99	248 - 260	216 - 228	29 - 35
5	67	101	254 - 266	222 - 234	30 - 36
6	68	104	261 - 273	229 - 241	31 - 37
7	70	106	267 - 279	235 - 247	32 - 38
8	71	109	274 - 286	242 - 254	32 - 38
9	73	111	280 - 292	248 - 260	33 - 39
10	74	114	287 - 299	255 - 267	34 - 40