

Men's - CLASSIC RANGE - ORIGINAL LAST - Fitting Chart D

Classic Range Original Lasts have a 9mm Inlay Allowance as Standard
 Footwear is supplied with 2 x 3mm Poron Inlays & 1 x Polyurethane Footbed as Standard
 (For 1/2 Sizes use measures half way between whole sizes)

Men's Slim Fitting

Regular Lasts						Adam Lasts (Narrow Heel)		Clive Lasts (Shallow Forepart)		Eddie Lasts (Narrow Heel & Shallow Forepart)			
Size	Heel Width (mm)	Joint Width (mm)	Joint Circ. (mm)	Instep Circ. (mm)	Toe Depth (mm)	Heel Width (mm)	Instep Circ. (mm)	Joint Circ. (mm)	Toe Depth (mm)	Heel Width (mm)	Joint Circ. (mm)	Instep Circ. (mm)	Toe Depth (mm)
1	54	77	198 - 216	192 - 210	29 - 38	50	185 - 203	191 - 209	25 - 34	50	191 - 209	185 - 203	25 - 34
2	55	80	205 - 223	199 - 217	29 - 38	51	192 - 210	198 - 216	25 - 34	51	198 - 216	192 - 210	25 - 34
3	57	82	211 - 229	205 - 223	30 - 39	53	199 - 217	205 - 223	26 - 35	53	205 - 223	199 - 217	26 - 35
4	58	85	218 - 236	212 - 230	31 - 40	54	205 - 223	211 - 229	27 - 36	54	211 - 229	205 - 223	27 - 36
5	60	87	224 - 242	218 - 236	31 - 40	56	212 - 230	218 - 236	27 - 36	56	218 - 236	212 - 230	27 - 36
6	61	90	231 - 249	225 - 243	32 - 41	57	218 - 236	224 - 242	28 - 37	57	224 - 242	218 - 236	28 - 37
7	63	92	237 - 255	231 - 249	33 - 42	59	225 - 243	231 - 249	29 - 38	59	231 - 249	225 - 243	29 - 38
8	64	95	244 - 262	238 - 256	33 - 42	60	231 - 249	237 - 255	29 - 38	60	237 - 255	231 - 249	29 - 38
9	66	97	250 - 268	244 - 262	34 - 43	62	238 - 256	244 - 262	30 - 39	62	244 - 262	238 - 256	30 - 39
10	67	100	257 - 275	251 - 269	35 - 44	63	244 - 262	250 - 268	31 - 40	63	250 - 268	244 - 262	31 - 40
11	69	102	263 - 281	257 - 275	35 - 44	65	251 - 269	257 - 275	31 - 40	65	257 - 275	251 - 269	31 - 40
12	70	105	270 - 288	264 - 282	36 - 45	66	257 - 275	263 - 281	32 - 41	66	263 - 281	257 - 275	32 - 41
13	72	107	276 - 294	270 - 288	37 - 46	68	264 - 282	270 - 288	33 - 42	68	270 - 288	264 - 282	33 - 42
14	73	110	283 - 301	277 - 295	37 - 46	69	270 - 288	276 - 294	33 - 42	69	276 - 294	270 - 288	33 - 42

Men's Standard Fitting

Regular Lasts						Adam Lasts (Narrow Heel)		Clive Lasts (Shallow Forepart)		Eddie Lasts (Narrow Heel & Shallow Forepart)			
Size	Heel Width (mm)	Joint Width (mm)	Joint Circ. (mm)	Instep Circ. (mm)	Toe Depth (mm)	Heel Width (mm)	Instep Circ. (mm)	Joint Circ. (mm)	Toe Depth (mm)	Heel Width (mm)	Joint Circ. (mm)	Instep Circ. (mm)	Toe Depth (mm)
1	58	83	210 - 228	200 - 218	29 - 38	54	193 - 211	203 - 221	25 - 34	54	203 - 221	193 - 211	25 - 34
2	59	86	217 - 235	207 - 225	29 - 38	55	200 - 218	210 - 228	25 - 34	55	210 - 228	200 - 218	25 - 34
3	61	88	223 - 241	213 - 231	30 - 39	57	207 - 225	217 - 235	26 - 35	57	217 - 235	207 - 225	26 - 35
4	62	91	230 - 248	220 - 238	31 - 40	58	213 - 231	223 - 241	27 - 36	58	223 - 241	213 - 231	27 - 36
5	64	93	236 - 254	226 - 244	31 - 40	60	220 - 238	230 - 248	27 - 36	60	230 - 248	220 - 238	27 - 36
6	65	96	243 - 261	233 - 251	32 - 41	61	226 - 244	236 - 254	28 - 37	61	236 - 254	226 - 244	28 - 37
7	67	98	249 - 267	239 - 257	33 - 42	63	233 - 251	243 - 261	29 - 38	63	243 - 261	233 - 251	29 - 38
8	68	101	256 - 274	246 - 264	33 - 42	64	239 - 257	249 - 267	29 - 38	64	249 - 267	239 - 257	29 - 38
9	70	103	262 - 280	252 - 270	34 - 43	66	246 - 264	256 - 274	30 - 39	66	256 - 274	246 - 264	30 - 39
10	71	106	269 - 287	259 - 277	35 - 44	67	252 - 270	262 - 280	31 - 40	67	262 - 280	252 - 270	31 - 40
11	73	108	275 - 293	265 - 283	35 - 44	69	259 - 277	269 - 287	31 - 40	69	269 - 287	259 - 277	31 - 40
12	74	111	282 - 300	272 - 290	36 - 45	70	265 - 283	275 - 293	32 - 41	70	275 - 293	265 - 283	32 - 41
13	76	113	288 - 306	278 - 296	37 - 46	72	272 - 290	282 - 300	33 - 42	72	282 - 300	272 - 290	33 - 42
14	77	116	295 - 313	285 - 303	37 - 46	73	278 - 296	288 - 306	33 - 42	73	288 - 306	278 - 296	33 - 42

Men's Wide Fitting

Regular Lasts						Adam Lasts (Narrow Heel)		Clive Lasts (Shallow Forepart)		Eddie Lasts (Narrow Heel & Shallow Forepart)			
Size	Heel Width (mm)	Joint Width (mm)	Joint Circ. (mm)	Instep Circ. (mm)	Toe Depth (mm)	Heel Width (mm)	Instep Circ. (mm)	Joint Circ. (mm)	Toe Depth (mm)	Heel Width (mm)	Joint Circ. (mm)	Instep Circ. (mm)	Toe Depth (mm)
1	62	89	222 - 240	208 - 226	29 - 38	58	201 - 219	215 - 232	25 - 34	58	215 - 232	201 - 219	25 - 34
2	63	92	229 - 247	215 - 233	29 - 38	59	208 - 226	222 - 240	25 - 34	59	222 - 240	208 - 226	25 - 34
3	65	94	235 - 253	221 - 239	30 - 39	61	215 - 233	229 - 247	26 - 35	61	229 - 247	215 - 233	26 - 35
4	66	97	242 - 260	228 - 246	31 - 40	62	221 - 239	235 - 253	27 - 36	62	235 - 253	221 - 239	27 - 36
5	68	99	248 - 266	234 - 252	31 - 40	64	228 - 246	242 - 260	27 - 36	64	242 - 260	228 - 246	27 - 36
6	69	102	255 - 273	241 - 259	32 - 41	65	234 - 252	248 - 266	28 - 37	65	248 - 266	234 - 252	28 - 37
7	71	104	261 - 279	247 - 265	33 - 42	67	241 - 259	255 - 273	29 - 38	67	255 - 273	241 - 259	29 - 38
8	72	107	268 - 286	254 - 272	33 - 42	68	247 - 265	261 - 279	29 - 38	68	261 - 279	247 - 265	29 - 38
9	74	109	274 - 292	260 - 278	34 - 43	70	254 - 272	268 - 286	30 - 39	70	268 - 286	254 - 272	30 - 39
10	75	112	281 - 299	267 - 285	35 - 44	71	260 - 278	274 - 292	31 - 40	71	274 - 292	260 - 278	31 - 40
11	77	114	287 - 305	273 - 291	35 - 44	73	267 - 285	281 - 299	31 - 40	73	281 - 299	267 - 285	31 - 40
12	78	117	294 - 312	280 - 298	36 - 45	74	273 - 291	287 - 305	32 - 41	74	287 - 305	273 - 291	32 - 41
13	80	119	300 - 318	286 - 304	37 - 46	76	280 - 298	294 - 312	33 - 42	76	294 - 312	280 - 298	33 - 42
14	81	122	307 - 325	293 - 311	37 - 46	77	286 - 304	300 - 318	33 - 42	77	300 - 318	286 - 304	33 - 42

Men's Deep Fitting

Regular Lasts					
Size	Heel Width (mm)	Joint Width (mm)	Joint Circ. (mm)	Instep Circ. (mm)	Toe Depth (mm)
1	58	83	222 - 240	200 - 218	35 - 44
2	59	86	229 - 247	207 - 225	35 - 44
3	61	88	235 - 253	213 - 231	36 - 45
4	62	91	242 - 260	220 - 238	37 - 46
5	64	93	248 - 266	226 - 244	37 - 46
6	65	96	255 - 273	233 - 251	38 - 47
7	67	98	261 - 279	239 - 257	39 - 48
8	68	101	268 - 286	246 - 264	39 - 48
9	70	103	274 - 292	252 - 270	40 - 49
10	71	106	281 - 299	259 - 277	41 - 50
11	73	108	287 - 305	265 - 283	41 - 50
12	74	111	294 - 312	272 - 290	42 - 51
13	76	113	300 - 318	278 - 296	43 - 52
14	77	116	307 - 325	285 - 303	43 - 52

Men's Wide & Deep Fitting

Regular Lasts					
Size	Heel Width (mm)	Joint Width (mm)	Joint Circ. (mm)	Instep Circ. (mm)	Toe Depth (mm)
1	62	89	234 - 252	208 - 226	35 - 44
2	63	92	241 - 259	215 - 233	35 - 44
3	65	94	247 - 265	221 - 239	36 - 45
4	66	97	254 - 272	228 - 246	37 - 46
5	68	99	260 - 278	234 - 252	37 - 46
6	69	102	267 - 285	241 - 259	38 - 47
7	71	104	273 - 291	247 - 265	39 - 48
8	72	107	280 - 298	254 - 272	39 - 48
9	74	109	286 - 304	260 - 278	40 - 49
10	75	112	293 - 311	267 - 285	41 - 50
11	77	114	299 - 317	273 - 291	41 - 50
12	78	117	306 - 324	280 - 298	42 - 51
13	80	119	312 - 330	286 - 304	43 - 52
14	81	122	319 - 337	293 - 311	43 - 52