

Paediatric - ALL LASTS - Fitting Chart F

Standard, Stability, Anti-Varus & Oversplint Lasts have a 3mm Inlay Allowance as Standard Footwear is supplied with 1 x 3mm Poron Inlay & 1 x Leather Covered Footbed as Standard

(Oversplint Lasts have 6mm extra width across midfoot & heels & 18mm extra instep circ. to accommodate AFO's)

Extra Slim Fitting							
Size (UK)	Size (EU)	Heel Width (mm)	Joint Width (mm)	Joint Circ. (mm)	Instep Circ. (mm)	Toe Depth (mm)	Leg Circ. (mm)
3	19	37	50	135	158	16	130
3 1/2		37	51	137	161	16	133
4	20	38	52	140	164	17	136
4 1/2		38	53	142	167	17	139
5	21	39	54	145	170	17	142
5 1/2	22	39	55	147	173	17	145
6	23	40	56	150	176	18	148
6 1/2		40	57	152	179	18	151
7	24	41	58	155	182	18	154
7 1/2	25	41	59	157	185	18	157
8		42	60	160	188	19	160
8 1/2	26	42	61	162	191	19	163
9	27	43	62	165	194	19	166
9 1/2		43	63	167	197	19	169
10	28	44	64	170	200	20	172
10 1/2	29	44	65	172	203	20	175
11		45	66	175	206	20	178
11 1/2	30	45	67	177	209	20	181
12	31	46	68	180	212	21	184
12 1/2		46	69	182	215	21	187
13	32	47	70	185	218	21	190
13 1/2		47	71	187	221	21	193
1	33	48	72	190	224	22	196
1 1/2	34	48	73	192	227	22	199
2		49	74	195	230	22	202
2 1/2	35	49	75	197	233	22	205
3	36	50	76	200	236	23	208
3 1/2		50	77	202	239	23	211
4	37	51	78	205	242	23	214
4 1/2		51	79	207	245	23	217
5	38	52	80	210	248	24	220

Slim Fitting							
Size (UK)	Size (EU)	Heel Width (mm)	Joint Width (mm)	Joint Circ. (mm)	Instep Circ. (mm)	Toe Depth (mm)	Leg Circ. (mm)
3	19	39	54	143	158	16	130
3 1/2		39	55	145	161	16	133
4	20	40	56	148	164	17	136
4 1/2		40	57	150	167	17	139
5	21	41	58	153	170	17	142
5 1/2	22	41	59	155	173	17	145
6	23	42	60	158	176	18	148
6 1/2		42	61	160	179	18	151
7	24	43	62	163	182	18	154
7 1/2	25	43	63	165	185	18	157
8		44	64	168	188	19	160
8 1/2	26	44	65	170	191	19	163
9	27	45	66	173	194	19	166
9 1/2		45	67	175	197	19	169
10	28	46	68	178	200	20	172
10 1/2	29	46	69	180	203	20	175
11		47	70	183	206	20	178
11 1/2	30	47	71	185	209	20	181
12	31	48	72	188	212	21	184
12 1/2		48	73	190	215	21	187
13	32	49	74	193	218	21	190
13 1/2		49	75	195	221	21	193
1	33	50	76	198	224	22	196
1 1/2	34	50	77	200	227	22	199
2		51	78	203	230	22	202
2 1/2	35	51	79	205	233	22	205
3	36	52	80	208	236	23	208
3 1/2		52	81	210	239	23	211
4	37	53	82	213	242	23	214
4 1/2		53	83	215	245	23	217
5	38	54	84	218	248	24	220

Standard Fitting							
Size (UK)	Size (EU)	Heel Width (mm)	Joint Width (mm)	Joint Circ. (mm)	Instep Circ. (mm)	Toe Depth (mm)	Leg Circ. (mm)
3	19	41	58	151	158	16	130
3 1/2		41	59	153	161	16	133
4	20	42	60	156	164	17	136
4 1/2		42	61	158	167	17	139
5	21	43	62	161	170	17	142
5 1/2	22	43	63	163	173	17	145
6	23	44	64	166	176	18	148
6 1/2		44	65	168	179	18	151
7	24	45	66	171	182	18	154
7 1/2	25	45	67	173	185	18	157
8		46	68	176	188	19	160
8 1/2	26	46	69	178	191	19	163
9	27	47	70	181	194	19	166
9 1/2		47	71	183	197	19	169
10	28	48	72	186	200	20	172
10 1/2	29	48	73	188	203	20	175
11		49	74	191	206	20	178
11 1/2	30	49	75	193	209	20	181
12	31	50	76	196	212	21	184
12 1/2		50	77	198	215	21	187
13	32	51	78	201	218	21	190
13 1/2		51	79	203	221	21	193
1	33	52	80	206	224	22	196
1 1/2	34	52	81	208	227	22	199
2		53	82	211	230	22	202
2 1/2	35	53	83	213	233	22	205
3	36	54	84	216	236	23	208
3 1/2		54	85	218	239	23	211
4	37	55	86	221	242	23	214
4 1/2		55	87	223	245	23	217
5	38	56	88	226	248	24	220

Wide Fitting							
Size (UK)	Size (EU)	Heel Width (mm)	Joint Width (mm)	Joint Circ. (mm)	Instep Circ. (mm)	Toe Depth (mm)	Leg Circ. (mm)
3	19	43	62	159	158	16	130
3 1/2		43	63	161	161	16	133
4	20	44	64	164	164	17	136
4 1/2		44	65	166	167	17	139
5	21	45	66	169	170	17	142
5 1/2	22	45	67	171	173	17	145
6	23	46	68	174	176	18	148
6 1/2		46	69	176	179	18	151
7	24	47	70	179	182	18	154
7 1/2	25	47	71	181	185	18	157
8		48	72	184	188	19	160
8 1/2	26	48	73	186	191	19	163
9	27	49	74	189	194	19	166
9 1/2		49	75	191	197	19	169
10	28	50	76	194	200	20	172
10 1/2	29	50	77	196	203	20	175
11		51	78	199	206	20	178
11 1/2	30	51	79	201	209	20	181
12	31	52	80	204	212	21	184
12 1/2		52	81	206	215	21	187
13	32	53	82	209	218	21	190
13 1/2		53	83	211	221	21	193
1	33	54	84	214	224	22	196
1 1/2	34	54	85	216	227	22	199
2		55	86	219	230	22	202
2 1/2	35	55	87	221	233	22	205
3	36	56	88	224	236	23	208
3 1/2		56	89	226	239	23	211
4	37	57	90	229	242	23	214
4 1/2		57	91	231	245	23	217
5	38	58	92	234	248	24	220