

Adjustable Cervical Collar



Indications

- Post surgery / post trauma
- Motor Neurone Disease and CVA
- Stable cervical fractures
- Cervical Spondylitis / Cervical Spondylosis
- O/A and R/A cervical spine
- Herniated cervical disc

Features

- X Ray and MRI compatible Velcro fastening cervical collar
- Versatile front and back components with a combination of low and high density poly ethylene shells allows stability and comfort
- Flexible tabs allow the collar to conform to the individual as it is tightened, distributing the pressure over the mandible and chest areas comfortably
- Occipital strap customises the collar at the rear to create an adjustable unique shelf to cushion and support the back of the head
- Breathable, cotton lined foam pads provide increased comfort and draws moisture away from the skin

Each Collar is supplied with a Replacement Liner

Size	Chin to shoulder (depth)	Code
Short	57mm or 2 1/4"	ACC-S
Regular	76mm or 3"	ACC-R
Tall	96mm or 3 1/4"	ACC-T
Extra Tall	114mm or 4 1/2"	ACC-XT



Chris Hanley & Partners
4 Kilvey Road, Brackmills, Northampton, NN4 7BQ
Tel: 01604 709999 Fax: 01604 700057