

EZ ROM Post Op Knee Brace



Features:

- Post-Op ROM brace with lightweight low profile hinge
- Ease of adjustment lets you make flexion and extension changes with a simple push of a button
- Compliance feature to help prevent tampering by adding user friendly security ties to lock the desired position
- Mouldable aluminium struts allow for custom contouring
- Aluminium strut length can also break down to accommodate sizing and rehab progression
- EZ ROM offers a full coverage lining for warmth and compression
- Universal fitting
- **NEW - Also now available as a COOL model - with 4 narrow fabric bands**

Indications:

- Post surgery e.g. TKR
- Post ligament repair
- Post fracture management
- Rehabilitation requiring ROM control
- Post trauma

Order Codes - Measure the circumference 6" above the mid-patella

FULL - Standard Length (24" 61cm)

Size	Circumference Range	Code
REG	Up to 27" (69cm)	EZR-S-R
X-LARGE	27-35" (69-89cm)	EZR-S-XL

FULL - Long Length (28" 71cm)

Size	Circumference Range	Code
REG	Up to 27" (69cm)	EZR-L-R
X-LARGE	27-35" (69-89cm)	EZR-L-XL

COOL - Standard Length (24" 61cm)

Size	Circumference Range	Code
REG	Up to 27" (69cm)	EZC-S-R
X-LARGE	27-35" (69-89cm)	EZC-S-XL

COOL - Long Length (28" 71cm)

Size	Circumference Range	Code
REG	Up to 27" (69cm)	EZC-L-R
X-LARGE	27-35" (69-89cm)	EZC-L-XL



Chris Hanley & Partners

4 Kilvey Road, Brackmills, Northampton, NN4 7BQ

Tel: 01604 709999 Fax: 01604 700057

e-mail: info@chaneco.co.uk www.chaneco.co.uk

