


## Finger Cots



### Features:

- Complete finger immobilisation
- Easy application
- Good finger tip protection
- Constructed of malleable aluminium
- With foam padding for patient comfort

### Indications:

- Post surgery
- Post trauma
- Finger fractures
- Tendon injuries
- Ligament injuries

### Order Codes

Size	Length (cm)	Code
X-Small	2.5	FC-XS
Small	3.8	FC-S
Medium	5.7	FC-M
Large	7.6	FC-L



**Chris Hanley & Partners**

**4 Kilvey Road, Brackmills, Northampton, NN4 7BQ**

**Tel: 01604 709999 Fax: 01604 700057**

**e-mail: [info@chaneco.co.uk](mailto:info@chaneco.co.uk) [www.chaneco.co.uk](http://www.chaneco.co.uk)**

