

## Heel Elevators



### Indications

- Shortening discrepancies of the lower limb
- Restricted ankle joint movement
- Achilles Tendonitis

### Features

- Firm Shore EVA Construction
- 3mm peel away layers
- Fabric cover

### Order Codes:

| Size         | Shoe Size (UK) | Shoe Size (Euro) | Code     |
|--------------|----------------|------------------|----------|
| Small/Medium | 3 - 8 1/2      | 36 - 42.5        | HE - S/M |
| Medium/Large | 9 - 12         | 43 - 47          | HE - M/L |



**Chris Hanley & Partners**

**4 Kilvey Road, Brackmills, Northampton, NN4 7BQ**

**Tel: 01604 709999 Fax: 01604 700057**

**e-mail: [info@chaneco.co.uk](mailto:info@chaneco.co.uk) [www.chaneco.co.uk](http://www.chaneco.co.uk)**

