

## Matrix Lite Medical Knee Brace



Silicone Pad



Knee Sleeve



Neoprene Wrap

### Indications

- ACL, PCL\*, MCL, LCL knee ligament instabilities
- Protection during low to medium impact activities
- \* Requires optional PCL Kit

### Features

- Strong, ultra lightweight (approximately 400g) frame constructed from heat adjustable carbon composite
- Shorter for greater mobility
- "ICR", low profile hinge mimics the natural movement of the knee
- Adjustable anterior tibial strap provides stability and adds support to the tibia
- Flexible cuffs, adjustable straps and condylar pads provide stability and customised fit
- Silicone pad in the lower cuff prevents migration and avoids tibial pressure
- Full ROM control via range of extension stops (included) and optional flexion stops (kit available separately)
- Shim kit for condylar pads and mini-screwdriver included for complete customised fitting
- Optional Knee Sleeve for added comfort, Neoprene Wrap for further customised fitting and replacement strap kit also available

### Measurements & Order Codes:

Size	Joint Width		Lite Medical		Knee Sleeve	Neoprene Wrap	PCL Kit	Flexion Stop Kit	Strap Kit
	Inches	Cm	Left	Right					
Small	3.5 - 4	8.9 - 10	MKB-LM-L-S	MKB-LM-R-S	MKB-KS-S	MKB-NSW	MKB-PCL	MKB-FSK	MKB-STR
Medium	4 - 4.5	10 - 11.4	MKB-LM-L-M	MKB-LM-R-M	MKB-KS-M				
Large	4.5 - 5	11.4 - 12.7	MKB-LM-L-L	MKB-LM-R-L	MKB-KS-L				
X Large	5 - 5.5	12.7 - 14	MKB-LM-L-XL	MKB-LM-R-XL	MKB-KS-XL				



Chris Hanley & Partners

4 Kilvey Road, Brackmills, Northampton, NN4 7BQ

Tel: 01604 709999 Fax: 01604 700057

e-mail: info@chaneco.co.uk www.chaneco.co.uk

