

Telescopic Elbow Brace



Indications

- Controlled elbow flexion/extension following surgery
- Trauma to elbow or muscles and ligaments of the arm
- Elbow arthroscopy
- Ulnar or radial collateral ligament reconstruction
- Stable or internally fixed fractures
- Medial / Lateral epicondylitis
- Severe sprains and strains

Features

- Telescopic struts and malleable aluminium cuffs lead to a more customised fit
- Single-sided lateral joint reduces bulk and allows for closer alignment to the arm
- Open, lightweight design, with ease of application and removal
- 3-dimensional padding and lycra lining offer superior comfort and compliance
- Optional shoulder strap has comfort pad at the neck and shoulder
- Flexion range can be quickly adjusted to 0, 10, 20, 30, 45, 60, 75, 90, 105 and 120°
- Extension can be adjusted to 0, 10, 20, 30, 45, 60, 75 and 90°
- Optional anti-supination/pronation arm bar kit is available

Model	Code
Left	TEB-L
Right	TEB-R
Arm Bar Kit	TEB-A

