

X-ROM 2



Indications

- Combined instabilities / injuries
- Condylar fractures
- Tibial plateau fractures
- Patella tendon repairs
- Osteochondral repairs
- Meniscal repairs
- High tibial osteotomy
- Total knees
- Knee sprains / strains

Features

- Easily adjustable telescoping hinge bars - from 44 - 69cm (17 - 27")
- ROM hinge allows 0°-120° flexion and -10°-+30° extension control, available in 10° increments
- One touch drop lock allows lock out at -10°, 0°, 10°, 20° and 30°
- Silicone coated breathable pads keep the patient cool and prevent migration
- Strap locking tabs to eliminate slippage
- Clip strap buckles for easy fastening
- Universal - one size fits all patients
- Lightweight (under 2lbs) brace increases patient compliance

Model	Code
Cool	XROM2C
Full	XROM2F

Price £84.95 + VAT

Free Delivery



Chris Hanley & Partners

4 Kilvey Road, Brackmills, Northampton, NN4 7BQ

Tel: 01604 709999 Fax: 01604 700057

e-mail: info@chaneco.co.uk www.chaneco.co.uk

