

# Smart STEPS Fitting Chart

Regular, Stability, Anti-Varus & Oversplint lasts all have a 3mm Inlay Allowance as standard.

Footwear is supplied with 1 x Leather covered footbed as standard

(Oversplint lasts have 8mm extra width across midfoot & heels & 20mm extra instep circ. to accommodate AFO's)

Infant Child Youth	Size (Uk)	Size (EU)	ALL FITTINGS			EXTRA SLIM FITTING			SLIM FITTING			STANDARD FITTING			WIDE FITTING		
			Instep Circ. (I)	Toe Depth (B)	Leg Circ. (H1)	Heel Width (D)	Joint Width (A)	Joint Circ. (J)	Heel Width (D)	Joint Width (A)	Joint Circ. (J)	Heel Width (D)	Joint Width (A)	Joint Circ. (J)	Heel Width (D)	Joint Width (A)	Joint Circ. (J)
			(mm)	(mm)	(mm)	(mm)	(mm)	(mm)	(mm)	(mm)	(mm)	(mm)	(mm)	(mm)	(mm)	(mm)	(mm)
I	2	18	155	15	127	36	49	133	38	53	141	40	57	149	42	61	157
I	3	19	158	18	130	37	50	135	39	54	143	41	58	151	43	62	159
I	3 ½		161	18	133	37	51	137	39	55	145	41	59	153	43	63	161
I	4	20	164	19	136	38	52	140	40	56	148	42	60	156	44	64	164
I	4 ½		167	19	139	38	53	142	40	57	150	42	61	158	44	65	166
I	5	21	170	19	142	39	54	145	41	58	153	43	62	161	45	66	169
I	5 ½	22	173	19	145	39	55	147	41	59	155	43	63	163	45	67	171
I	6	23	176	20	148	40	56	150	42	60	158	44	64	166	46	68	174
I	6 ½		179	20	151	40	57	152	42	61	160	44	65	168	46	69	176
I	7	24	182	20	154	41	58	155	43	62	163	45	66	171	47	70	179
I	7 ½	25	185	20	157	41	59	157	43	63	165	45	67	173	47	71	181
I	8		188	21	160	42	60	160	44	64	168	46	68	176	48	72	184
I	8 ½	26	191	21	163	42	61	162	44	65	170	46	69	178	48	73	186
I	9	27	194	21	166	43	62	165	45	66	173	47	70	181	49	74	189
I	9 ½		197	21	169	43	63	167	45	67	175	47	71	183	49	75	191
I	10	28	200	22	172	44	64	170	46	68	178	48	72	186	50	76	194
C	10 ½	29	203	22	175	44	65	172	46	69	180	48	73	188	50	77	196
C	11		206	22	178	45	66	175	47	70	183	49	74	191	51	78	199
C	11 ½	30	209	22	181	45	67	177	47	71	185	49	75	193	51	79	201
C	12	31	212	23	184	46	68	180	48	72	188	50	76	196	52	80	204
C	12 ½		215	23	187	46	69	182	48	73	190	50	77	198	52	81	206
C	13	32	218	23	190	47	70	185	49	74	193	51	78	201	53	82	209
C	13 ½		221	23	193	47	71	187	49	75	195	51	79	203	53	83	211
C	1	33	224	24	196	48	72	190	50	76	198	52	80	206	54	84	214
C	1 ½	34	227	24	199	48	73	192	50	77	200	52	81	208	54	85	216
C	2		230	24	202	49	74	195	51	78	203	53	82	211	55	86	219
C	2 ½	35	233	24	205	49	75	197	51	79	205	53	83	213	55	87	221
C	3	36	236	25	208	50	76	200	52	80	208	54	84	216	56	88	224
C	3 ½		239	25	211	50	77	202	52	81	210	54	85	218	56	89	226
C	4	37	242	25	214	51	78	205	53	82	213	55	86	221	57	90	229
C	4 ½		245	25	217	51	79	207	53	83	215	55	87	223	57	91	231
C	5	38	248	26	220	52	80	210	54	84	218	56	88	226	58	92	234
Y	5 ½	39	251	26	223	52	81	212	54	85	220	56	89	228	58	93	236
Y	6		254	26	226	53	82	215	55	86	223	57	90	231	59	94	239
Y	6 ½	40	257	26	229	53	83	217	55	87	225	57	91	233	59	95	241
Y	7	41	260	26	232	54	84	220	56	88	228	58	92	236	60	96	244
Y	7 ½		263	27	235	54	85	222	56	89	230	58	93	238	60	97	246
Y	8	42	266	27	238	55	86	225	57	90	233	59	94	241	61	98	249
Y	8 ½		269	27	241	55	87	227	57	91	235	59	95	243	61	99	251
Y	9	43	272	27	244	56	88	230	58	92	238	60	96	246	62	100	254
Y	9 ½	44	275	27	247	56	89	232	58	93	240	60	97	248	62	101	256