



In addition to the standard inlays supplied in the footwear, we are also able to offer a range of alternative inlays manufactured from heat mouldable EVA with a microfibre top cover to fit our footwear perfectly. All inlays are available in all sizes and width fittings with heel cups up to a height of 25mm (please specify when ordering).



## Neuro Inlays

- Supports valgus and varus arches
- Neutral heel cup
- Stabilises the metatarsal joints with a moulded metatarsal dome and incorporated dropped hallux
- Wedges, posts and raises can be added for further control



## Neutral Inlays

- Supports valgus and varus arches
- Deep neutral heel cup
- Helps to control overpronation by reducing strain and pressure on arches
- Wedges, posts and raises can be added for further control



## Anti Pronation Inlays

- Supports valgus and varus arches
- Incorporates a kirby skive and neutral extrinsic post
- Corrects and controls rearfoot
- Wedges, posts and raises can be added for further control

Our standard Arch and Filled Arch Support Inlays both offer extra support to maintain a good biomechanical foot position to help prevent fallen arches, plantar fasciitis and further damage.



## Arch Support Inlays

- Available in full and  $\frac{3}{4}$  lengths
- Supports valgus and varus arches
- Controls moderate overpronation
- Valgus & varus arches can be made up to a height of 25mm
- Wedges, posts and raises can be added for further control



## Filled Arch Support Inlays

- Available in full and  $\frac{3}{4}$  lengths
- Supports valgus and varus arches
- Controls higher level overpronation with filled arches
- Valgus & varus arches can be made up to a height of 25mm
- Wedges, posts and raises can be added for further control