



Smart  
STEPS

# Lasts & Soles




## Regular Last




All styles will be made on our Regular last with standard height long medial and lateral stiffeners unless specified otherwise

## Stability Option




The stability option is available on all styles made on our Regular or Anti-Varus lasts and offers a higher and longer stiffener for increased support

## Anti-Varus Last



All styles can be made on our Anti-Varus last incorporating 8° forefoot abduction for the treatment of adducted or club feet

## Oversplint Last



All styles can be made on our Oversplint last for the accommodation of AFO's providing a straight heel counter, a wider rearfoot and increased instep & ankle circumference

## Sole Units

All stitchdown footwear will be supplied with a hard wearing rubber commando sole in either black or caramel. We also offer a lightweight EVA sole should weight of the footwear be paramount above wearability



Rubber Commando Sole



Lightweight EVA Sole



All our trainer styles have a cemented construction with either a Side Stitched sole or Trainer Sole depending on the style in either black or white. Both sole units have a rubber bottom which could be changed to lightweight EVA to reduce weight



Side Stitched Sole



Trainer Sole