



FITTING INSTRUCTIONS

ROM Elbow Brace

In stock and available for next working day delivery



Features:

- Aluminium ROM hinge allows flexion and extension adjustments
- Telescoping struts allows length to be adjusted from 13" to 17" for a customised fit
- Adjustable hook and loop sling straps provide additional immobilisation
- Non-slip, breathable foam liner
- Open, lightweight design with ease of application and removal

Indications:

- Post-operative / Post-trauma support
- Elbow arthroscopy
- Ulnar or radial collateral ligament reconstruction
- Stable or internally fixed fractures
- Medial / Lateral epicondylitis
- Severe strains and sprains

Side

Code

NHSSC

Left

REB-L

GRZ964

Right

REB-R

GRZ965