



NINEWELLS BOOTS

Ninewells Boots

Fully adjustable unisex boot, opens to toe for maximum volume



Features:

- Opens to toe for easy access, unlined
- 3 hook and loop straps for complete adjustability and maximum volume
- Wide and Deep fitting
- Through Wedge Rocker sole
- Available in with black brushed nylon or black neoprene, in sizes 1 - 14
- Sold in pairs

Indications:

- To accommodate bandaged feet, oedema, deformities and prostheses
- Diabetic foot problems

Also available as a made-to-order item in half sizes and for any adaptations required

Part of the Chaneco Multi-Footwear Range:



Multi Boot



Multi Shoe



Multi Sandal



Ultra Boot



Chaneco



4 Kilvey Road
Brackmills
Northampton
NN4 7BQ

T: 01604 709999

E: info@chaneco.co.uk

Ninewells Boots

Order Codes

UK Shoe Size	Brushed Nylon Chaneco Code	Brushed Nylon NHSSC Code	Neoprene Chaneco Code	Neoprene NHSSC Code
1	NWB-1	ERP5336	NWBN-1	N/A
2	NWB-2	ERP5342	NWBN-2	ERP5393
3	NWB-3	ERP5343	NWBN-3	ERP5394
4	NWB-4	ERP5344	NWBN-4	ERP5395
5	NWB-5	ERP5345	NWBN-5	ERP5396
6	NWB-6	ERP5383	NWBN-6	ERP5397
7	NWB-7	ERP5384	NWBN-7	ERP5398
8	NWB-8	ERP5385	NWBN-8	ERP5399
9	NWB-9	ERP5386	NWBN-9	ERP5400
10	NWB-10	ERP5337	NWBN-10	ERP5387
11	NWB-11	ERP5338	NWBN-11	ERP5388
12	NWB-12	ERP5339	NWBN-12	ERP5389
13	NWB-13	ERP5340	NWBN-13	ERP5390
14	NWB-14	ERP5341	NWBN-14	ERP5391