

NEW Products....

We are excited to let you know about these latest additions to the Chaneco orthotics range from our supplier OCSI, the leading provider of contracture management orthotic devices in the US. All now in stock and available to order today.

For more information and to discuss pricing, please contact your Chaneco Business Manager. Scan the QR codes below for contact details or call our Customer Services Team on 01604 709999.



Elite Rehabilitator AFO

A carbon fibre AFO with a unique flexible footplate to facilitate hip and knee flexion for gait rehabilitation and relearning.



Size	Shoe Size		Code	
	UK	Euro	Left	Right
X Small	<1.5	<33	AFO-ER-XS-L	AFO-ER-XS-R
Small	2 - 4	34 - 36.5	AFO-ER-S-L	AFO-ER-S-R
Medium	4.5 - 7	37 - 41	AFO-ER-M-L	AFO-ER-M-R
Large	7.5 - 10	42 - 44	AFO-ER-L-L	AFO-ER-L-R
X Large	10.5+	45+	AFO-ER-XL-L	AFO-ER-XL-R

Elite Rehabilitator Pro AFO

The Pro version of the Elite Rehabilitator has an enhanced anterior shell and lateral strut, providing more significant support.



Size	Shoe Size		Code	
	UK	Euro	Left	Right
X Small	<1.5	<33	AFO-ERP-XS-L	AFO-ERP-XS-R
Small	2 - 4	34 - 36.5	AFO-ERP-S-L	AFO-ERP-S-R
Medium	4.5 - 7	37 - 41	AFO-ERP-M-L	AFO-ERP-M-R
Large	7.5 - 10	42 - 44	AFO-ERP-L-L	AFO-ERP-L-R
X Large	10.5	45+	AFO-ERP-XL-L	AFO-ERP-XL-R

Elite Rehabilitator Spiral AFO

The spiral lateral strut with posterior calf cuff provides superior gait assistance for more active patients needing dynamic dorsi-assist to correct foot drop.



Size	Shoe Size		Code	
	UK	Euro	Left	Right
X Small	<1.5	<33	AFO-ERS-XS-L	AFO-ERS-XS-R
Small	2 - 4	34 - 36.5	AFO-ERS-S-L	AFO-ERS-S-R
Medium	4.5 - 7	37 - 41	AFO-ERS-M-L	AFO-ERS-M-R
Large	7.5 - 10	42 - 44	AFO-ERS-L-L	AFO-ERS-L-R
X Large	10.5+	45+	AFO-ERS-XL-L	AFO-ERS-XL-R

South



Stuart Mills

Midlands



Ian Clifton

North



Tracey Brown


Scotland / N. Ireland



Darren Loughrey


4 Kilvey Road, Brackmills
Northampton, NN4 7BQ
T: 01604 709999
E: info@chaneco.co.uk

DynaPro Torticollis Orthosis



Clinically indicated to provide low load prolonged stretch to the neck to reverse Torticollis. Improves range of motion and functional alignment of the head and neck.

DynaPro Finger Flex



Resting hand design with a heat mouldable thermoplastic insert to provide therapeutic neuro-inhibiting stretch for the treatment of contractures.

Size	Forehead to Shoulder		Code	
	Inches	CM	Left	Right
Regular	6 - 7	15 - 18	DPTO-R-L	DPTO-R-R
Large	7 - 9	18 - 23	DPTO-L-L	DPTO-L-R

Size	Length CM	MP Joint Width		Code	
		Inches	CM	Left	Right
Small	28.5	2½ - 3	6 - 7.5	DPFF-S-L	DPFF-S-R
Medium	31	3 - 3½	7.5 - 9	DPFF-M-L	DPFF-M-R
Large	33.5	3½ - 4	9 - 10	DPFF-L-L	DPFF-L-R

DynaPro Flex Elbow



Ideal "static-dynamic" therapy to reverse mild to moderate tissue shortening contractures of the elbow, particularly suitable for patients with abnormal tone and spasticity. Modifiable over time as rehabilitation progresses.

OrthoPro HyperEx Knee

Size	Thigh Circumference		Code	
	Inches	CM	Left	Right
Small	9 - 13	23 - 33	KOP-S-L	KOP-S-R
Medium	14 - 19	35.5 - 48	KOP-M-L	KOP-M-R
Large	20 - 25	51 - 63.5	KOP-L-L	KOP-L-R
X Large	26 - 31	66 - 79	KOP-XL-L	KOP-XL-R

Effective control of genu recurvatum via an adjustable pneumatic air bladder behind the knee joint, improving gait by encouraging knee flexion. Easy to wear whilst walking and sitting.



Size	Lower Arm Length - CM	Upper Arm Circ. Inches	Upper Arm Circ. CM	Code
Small	12.5 - 20	8 - 10	20 - 25.5	DPE-S
Medium	19 - 25.5	9½ - 11½	24 - 29	DPE-M
Large	23 - 30.5	11 - 14	28 - 35.5	DPE-L

Scan here for a digital copy of our Orthotics catalogue:



ORTHOTICS

4 Kilvey Road, Brackmills
Northampton, NN4 7BQ
T: 01604 709999
E: info@chaneco.co.uk