

Silicone Gel Toe Cup
- Fitting Guidelines and Care Instructions



Features:

- Pre-formed silicone gel toe cup
- Covers and protects the end of the toe
- Sold singularly

Indications:

- All round protection of the toe

Fitting Instructions:

- Use as per clinician's guidelines

Care & Maintenance:

- Hand wash in warm water with mild soap
- Allow to dry naturally away from direct heat sources
- Do not use bleach

Size Table

Size	Code
Small/Medium	STC-S/M
Large/X-Large	STC-L/XL

CHANECO LTD
 4 Kilvey Road, Brackmills, Northampton, NN4 7BQ
 T: 01604 709999 F: 01604 700057
 www.chaneco.co.uk



STC IFU - Revision 2 - Jan 2024



Chaneco Limited is registered in England under Company Number 11511502. Registered office as above.