

Silicone Toe Crest Cushion
- Fitting Guidelines and Care Instructions



- Features:**
- Silicone toe prop which fit snugly between the ball of the foot and the toes
 - Relieves discomfort by raising and realigning the toes
 - Eases pressure on metatarsal area of bunching toes
 - Sold in pairs

- Indications:**
- Hammer, claw or mallet toes
 - General toe deformities

- Fitting Instructions:**
- Use as per clinician's guidelines

- Care & Maintenance:**
- Hand wash in warm water with mild soap
 - Allow to dry naturally away from direct heat sources

Size Table

Size	Length (mm)	Code
Small	63	STCC-S
Medium	70	STCC-M
Large	77	STCC-L

CHANECO LTD
4 Kilvey Road, Brackmills, Northampton, NN4 7BQ
T: 01604 709999 F: 01604 700057
www.chaneco.co.uk

Chaneco Limited is registered in England under Company Number 11511502. Registered office as above.



STCC IFU - Revision 2 - Jan 2024

