

Silicone Gel Toe/Corn Shield
- Fitting Guidelines and Care Instructions



Features:

- 2mm thick gel disc protects toes from pressure and friction
- Outside covered to reduce sheer
- Sold singularly

Indications:

- Pressure areas on clawing toes
- Pain relief from corns, calluses or blisters
- Post-trauma/surgery

Fitting Instructions:

- Use as per clinician's guidelines

Care & Maintenance:

- Hand wash in warm water with mild soap
- Allow to dry naturally away from direct heat sources
- Do not use bleach

Size Table

Size	Width		Code
	Inches	CM	
Small	1	2.5	STCS-S
Medium	1½	3	STCS-M
Large	1¾	3.5	STCS-L

CHANECO LTD
 4 Kilvey Road, Brackmills, Northampton, NN4 7BQ
 T: 01604 709999 F: 01604 700057
 www.chaneco.co.uk



STCS IFU - Revision 2 - Jan 2024



Chaneco Limited is registered in England under Company Number 11511502. Registered office as above.