

Silicone Toe Spreader
- Fitting Guidelines and Care Instructions



- Features:**
- Silicone corrective shield between the toes
 - Reduces friction by keeping toes apart
 - Slimmer in the middle for greater comfort with broader top and bottom to maintain spread
 - Sold singularly


- Indications:**
- General toe deformities
 - Post-surgery to maintain alignment

- Fitting Instructions:**
- Use as per clinician's guidelines

- Care & Maintenance:**
- Hand wash in warm water with mild soap
 - Allow to dry naturally away from direct heat sources

Size Table

| Size | Height (mm) | Code |
|--------|-------------|--------|
| Small | 20 | STSP-S |
| Medium | 25 | STSP-M |
| Large | 30 | STSP-L |

 **CHANECO LTD**
4 Kilvey Road, Brackmills, Northampton, NN4 7BQ
T: 01604 709999 F: 01604 700057
www.chaneco.co.uk
Chaneco Limited is registered in England under Company Number 11511502. Registered office as above.



STSP IFU - Revision 2 - Jan 2024

