

Silicone Gel Tubular Toe Bandage
- Fitting Guidelines and Care Instructions



Features:

- Elasticated tube internally coated throughout with silicone gel
- All round protection from pressure and friction
- Tube length 145mm (5 7/8") - cut to desired length
- Outside covered to reduce sheer

Indications:

- Pressure areas on clawing toes
- Pain relief from corns, calluses or blisters
- Post-trauma/surgery

Fitting Instructions:

- Use as per clinician's guidelines

Care & Maintenance:

- Hand wash in warm water with mild soap
- Allow to dry naturally away from direct heat sources
- Do not use bleach

Size Table

| Size | Width | | Code |
|--------|--------|-----|--------|
| | Inches | CM | |
| Small | 1 | 2.5 | STTB-S |
| Medium | 1 1/4 | 3 | STTB-M |
| Large | 1 1/2 | 3.5 | STTB-L |

STTB IFU - Revision 2 - Jan 2024

CHANECO LTD
 4 Kilvey Road, Brackmills, Northampton, NN4 7BQ
 T: 01604 709999 F: 01604 700057
 www.chaneco.co.uk
 Chaneco Limited is registered in England under Company Number

